

Disposable protective masks - 10 Box of 50 48.00 €Ex VAT Covid protection material







FEATURES

Disposable 3-ply surgical protective mask

Brand Base Fitting Ref ID Delivery schedule Category Type Securite shopping

No fixing mas-de-pro-3071 3071 1-3 working days Covid protection material Security system

DETAIL

Pack of 10 boxes of 50 disposable surgical masks.

Surgical masks type II - ADULT

3 ply - non sterile and single use Bacterial filtration efficiency> 98% (type 2).

Filter droplets, pollen, dust.

Composed of 3 layers (initiated support / polypropylene fiber filter).

Covers nose and mouth.

DESCRIPTION

Ear fixation by elastic.

Packaging and legal notices in FRENCH).

Blue non-woven mask, one size with elastic for better protection. To prevent the spread of Covid-19, colds and flu, it protects both the wereine of Energy against pathogenic germs and those around him from any microbial dissemination.

To ensure maximum protection, the mask must be replaced regularly during the day. The recommended wearing time is 4 hours maximum.

Disposable protective masks - 10 Box of 50 Disposable protective masks - 10 boxes or 50 type 2 adult surgical masks - Complex with standard EN14683 / 2019



FEATURES

Brand Base Fitting Ref ID Delivery schedule Category Type Securite shopping

No fixing mas-de-pro-3071 3071 1-3 working days Covid protection material Security system

DETAIL

Pack of 10 boxes of 50 disposable surgical masks.

Surgical masks type II - ADULT

3 ply - non sterile and single use Bacterial filtration efficiency> 98% (type 2).

Filter droplets, pollen, dust.

Composed of 3 layers (initiated support / polypropylene fiber filter).

.....

Covers nose and mouth.

Disposable protective masks - Dimensions of the surgical mask: 17 cm x 9.5 cm. Covid protection material - Photo n° 2 Packaging and logal notices in EPEN(CH)



Descriptif de nos masques chirurgicaux :

- non stérile et à usage unique

- 3 plis
- ≥98 % d'efficacité de filtration bactérienne
- Filtre gouttelettes, pollen, poussières
- Composé de 3 épaisseurs
- couvre le nez et la bouche
- Dimensions du masque chirurgical : 17 cm x 9,5 cm
- Fixation auriculaire par élastique
- Matière : support en intissé, filtre en fibre de Polypropylène
- -Durée de conservation : 2 ans
- -Disponible en boite de 50 masques



Nos masques sont certifiés conforme à la norme européenne en vigueur EN14683 : 2019

FEATURES		DETAIL		
Brand	Securite shopping	Pack of 10 boxes of 50 disposable surgical masks.		
Base		Surgical masks type II - ADULT		
Fitting	No fixing	5 51		
Ref	mas-de-pro-3071	3 ply - non sterile and single use		
ID	3071	Bacterial filtration efficiency> 98% (type 2).		
Delivery schedule Category	1-3 working days Covid protection material	Filter droplets, pollen, dust.		
Туре	Security system	Composed of 3 layers (initiated support / polypropylene fiber filter).		
		Covers nose and mouth.		

Disposable protective masks - Dimensions of the surgical mask: 17 cm x 9.5 cm. Covid protection material - Photo n° 3 Fiche technique des masques chirurgicaux

Pour une protection optimale, le masque chirurgical doit être posé avec attention sur votre visage. Les étapes suivantes vous permettrons d'être protéger efficacement :

O- Se laver les mains et placer le haut du masque sur le nez et les bandes élastiques autour des oreilles, face bleue à l'extérieur

2 - Ajuster la baguette du bord haut du masque sur le nez pour adapter la forme à votre morphologie, puis déployer le bas du masque sous votre menton

3 - Procéder aux derniers ajustements pour assurer une couverture efficace du nez et de la bouche



6 - Lavez-vous les mains après avoir rétiré votre masque, et jetez-le dans une poubelle fermée après usage

Mise en garde :

- Ce produit est à usage unique, il est jetable et ne peut pas être lavé et réutilisé.
- Durée d'utilisation recommandée : 4 heures
- Ne pas utiliser en salle stérile ou en salle d'opération
- Ce produit ne doit pas être utilisé après sa date limite de conservation (voir sur le dessous de la boite ou si l'emballage est endommagé)
- Le produit doit être stocké au sec et à l'abri de température extrêmes





F	EΑ	UF	RES

BrandSecurite shoppingBaseNo fixingFittingNo fixingRefmas-de-pro-3071ID3071Delivery schedule1-3 working daysCategoryCovid protection materialTypeSecurity system

DETAIL

Pack of 10 boxes of 50 disposable surgical masks.

Surgical masks type II - ADULT

3 ply - non sterile and single use Bacterial filtration efficiency> 98% (type 2).

Filter droplets, pollen, dust.

Composed of 3 layers (initiated support / polypropylene fiber filter).

Covers nose and mouth.

Disposable protective masks - Dimensions of the sprotectionmask: 17 cm x 9.5 cm. Covid protection material - Photo n 4



FEATURES

Brand Base Fitting Ref ID Delivery schedule Category Type Securite shopping No fixing mas-de-pro-3071 3071 1-3 working days Covid protection material Security system

DETAIL

Pack of 10 boxes of 50 disposable surgical masks.

Surgical masks type II - ADULT

3 ply - non sterile and single use Bacterial filtration efficiency> 98% (type 2).

Filter droplets, pollen, dust.

Composed of 3 layers (initiated support / polypropylene fiber filter).

Covers nose and mouth.

Disposable protective masks - Dimensions of the surgical mask: 17 cm x 9.5 cm. **Covid protection material** - Photo n 5 Packaging and logal patient in EPENICH



FEATURES

Brand Base Fitting Ref ID Delivery schedule Category Type

Securite shopping No fixing

mas-de-pro-3071 3071 1-3 working days Covid protection material Security system

DETAIL

Pack of 10 boxes of 50 disposable surgical masks.

Surgical masks type II - ADULT

3 ply - non sterile and single use Bacterial filtration efficiency> 98% (type 2).

Filter droplets, pollen, dust.

Composed of 3 layers (initiated support / polypropylene fiber filter).

Covers nose and mouth.

Dimensions of the surgical mask: 17 cm x 9.5 cm.

Ear fixation by elastic.

Packaging and legal notices in FRENCH).

CE certified - EN14683: 2019